

SMALL PACK SILICA GEL PORTFOLIO

Material Type: Granular

CAS No.: 7631-86-9

REACH No.: 01-2119379799-16-0299

Cat. No.	Mfr. No.	Description	Pore Size	Particle Size	Derivatisation	Pack Size
30048966	301989	ZEOprep 60 / 25-40µm	60Å	25-40µm	2-4%	1kg
30048948	301988	ZEOprep 60 / 25-40µm	60Å	25-40µm	2-4%	5kg
30049033	300231	ZEOprep 60 / 40-63µm	60Å	40-63µm	2-4%	1kg
30048995	302097	ZEOprep 60 / 40-63µm	60Å	40-63µm	2-4%	5kg
30048958	300227	ZEOprep 60 / 60-200µm	60Å	60-200µm	2-4%	1kg
30048987	302098	ZEOprep 60 / 60-200µm	60Å	60-200µm	2-4%	5kg
30048972	302099	ZEOprep 60 ACT / 35-70µm	60Å	35-70µm	<1%	1kg
30048960	302100	ZEOprep 60 ACT / 35-70µm	60Å	35-70µm	<1%	5kg
30049010	300214	ZEOprep 60 HYD / 40-63µm	60Å	40-63µm	4-6%	1kg
30048982	302103	ZEOprep 60 HYD / 40-63µm	60Å	40-63µm	4-6%	5kg
30048990	300211	ZEOprep 60 HYD / 60-200µm	60Å	60-200µm	4-6%	1kg
30048953	302104	ZEOprep 60 HYD / 60-200µm	60Å	60-200µm	4-6%	5kg
30049013	302105	ZEOprep 90 / 40-63µm	90Å	40-63µm	2-4%	1kg
30048952	302106	ZEOprep 90 / 40-63µm	90Å	40-63µm	2-4%	5kg
30049009	302107	ZEOprep 90 / 60-200µm	90Å	60-200µm	2-4%	1kg
30048976	302108	ZEOprep 90 / 60-200µm	90Å	60-200µm	2-4%	5kg

SMALL PACK SILICA GEL PORTFOLIO

Material Type: Functionised Granular

CAS No.: 7631-86-9

REACH No.: 01-2119379799-16-0299

Cat. No.	Mfr. No.	Description	Pore Size	Particle Size	Derivatisation	Pack Size
30048957	300190	ZEOprep 60 C ₁₈ / 40-63µm	60Å	40-63µm	C ₁₈	100g
30048967	302074	ZEOprep 60 C ₁₈ / 40-63µm	60Å	40-63µm	C ₁₈	500g
30049024	300191	ZEOprep 60 C ₁₈ / 40-63µm	60Å	40-63µm	C ₁₈	1kg
30049031	302089	ZEOprep 60 C ₁₈ / 40-63µm	60Å	40-63µm	C ₁₈	5kg
30048971	300175	ZEOprep 90 C ₁₈ / 40-63µm	90Å	40-63µm	C ₁₈	100g
30049006	302075	ZEOprep 90 C ₁₈ / 40-63µm	90Å	40-63µm	C ₁₈	500g
30049017	300176	ZEOprep 90 C ₁₈ / 40-63µm	90Å	40-63µm	C ₁₈	1kg
30048950	302090	ZEOprep 90 C ₁₈ / 40-63µm	90Å	40-63µm	C ₁₈	5kg
30048965	300186	ZEOprep 60 DIOL / 40-63µm	60Å	40-63µm	Diol	100g
30048949	302076	ZEOprep 60 DIOL / 40-63µm	60Å	40-63µm	Diol	500g
30049026	300187	ZEOprep 60 DIOL / 40-63µm	60Å	40-63µm	Diol	1kg
30049002	302091	ZEOprep 60 DIOL / 40-63µm	60Å	40-63µm	Diol	5kg
30048962	300184	ZEOprep 60 NH ₂ / 40-63µm	60Å	40-63µm	NH ₂	100g
30048955	302080	ZEOprep 60 NH ₂ / 40-63µm	60Å	40-63µm	NH ₂	500g
30048997	300185	ZEOprep 60 NH ₂ / 40-63µm	60Å	40-63µm	NH ₂	1kg
30048973	302092	ZEOprep 60 NH ₂ / 40-63µm	60Å	40-63µm	NH ₂	5kg
30048991	300182	ZEOprep 60 NH ₂ P / 40-63µm	60Å	40-63µm	NH ₂ P	100g
30048984	302081	ZEOprep 60 NH ₂ P / 40-63µm	60Å	40-63µm	NH ₂ P	500g
30049020	300183	ZEOprep 60 NH ₂ P / 40-63µm	60Å	40-63µm	NH ₂ P	1kg
30048963	302093	ZEOprep 60 NH ₂ P / 40-63µm	60Å	40-63µm	NH ₂ P	5kg
30049021	300166	ZEOprep 60 SAX / 40-63µm	60Å	40-63µm	SAX	100g
30048988	302082	ZEOprep 60 SAX / 40-63µm	60Å	40-63µm	SAX	500g
30049015	300167	ZEOprep 60 SAX / 40-63µm	60Å	40-63µm	SAX	1kg
30049007	302094	ZEOprep 60 SAX / 40-63µm	60Å	40-63µm	SAX	5kg
30048983	300177	ZEOprep 60 SH / 40-63µm	60Å	40-63µm	SH	100g
30049005	302083	ZEOprep 60 SH / 40-63µm	60Å	40-63µm	SH	500g
30049030	300178	ZEOprep 60 SH / 40-63µm	60Å	40-63µm	SH	1kg
30049022	302095	ZEOprep 60 SH / 40-63µm	60Å	40-63µm	SH	5kg

SMALL PACK SILICA GEL PORTFOLIO

Material Type: Spheroidal

CAS No.: 7631-86-9

REACH No.: 01-2119379799-16-0299

Cat. No.	Mfr. No.	Description	Pore Size	Particle Size	Derivatisation	Pack Size
30049029	300257	ZEObeads 50 C ₁₈ / 20-30µm	50Å	20-30µm	C ₁₈	100g
30048974	302070	ZEObeads 50 C ₁₈ / 20-30µm	50Å	20-30µm	C ₁₈	500g
30048961	300258	ZEObeads 50 C ₁₈ / 20-30µm	50Å	20-30µm	C ₁₈	1kg
30048999	302085	ZEObeads 50 C ₁₈ / 20-30µm	50Å	20-30µm	C ₁₈	5kg
30048978	300255	ZEObeads 50 C ₈ / 20-30µm	50Å	20-30µm	C ₈	100g
30048954	302071	ZEObeads 50 C ₈ / 20-30µm	50Å	20-30µm	C ₈	500g
30048959	300256	ZEObeads 50 C ₈ / 20-30µm	50Å	20-30µm	C ₈	1kg
30049014	302086	ZEObeads 50 C ₈ / 20-30µm	50Å	20-30µm	C ₈	5kg
30049025	300251	ZEObeads 50 NH ₂ / 20-30µm	50Å	20-30µm	NH ₂	100g
30048992	302072	ZEObeads 50 NH ₂ / 20-30µm	50Å	20-30µm	NH ₂	500g
30049000	300252	ZEObeads 50 NH ₂ / 20-30µm	50Å	20-30µm	NH ₂	1kg
30048985	302087	ZEObeads 50 NH ₂ / 20-30µm	50Å	20-30µm	NH ₂	5kg
30048989	300253	ZEObeads 50 DIOL / 20-30µm	50Å	20-30µm	Diol	100g
30049003	302073	ZEObeads 50 DIOL / 20-30µm	50Å	20-30µm	Diol	500g
30048968	300254	ZEObeads 50 DIOL / 20-30µm	50Å	20-30µm	Diol	1kg
30048956	302088	ZEObeads 50 DIOL / 20-30µm	50Å	20-30µm	Diol	5kg

SMALL PACK SILICA GEL PORTFOLIO

Material Type: Spherical

CAS No.: 7631-86-9

REACH No.: 01-2119379799-16-0299

Cat. No.	Mfr. No.	Description	Pore Size	Particle Size	Derivatisation	Pack Size
30048964	300242	ZEOsphere 100 C ₁₈ / 10µm	100Å	10µm	C ₁₈	100g
30049019	302054	ZEOsphere 100 C ₁₈ / 10µm	100Å	10µm	C ₁₈	500g
30049027	300933	ZEOsphere 100 C ₁₈ / 10µm	100Å	10µm	C ₁₈	1kg
30049012	300270	ZEOsphere 120 C ₁₈ / 10µm	120Å	10µm	C ₁₈	100g
30048996	302055	ZEOsphere 120 C ₁₈ / 10µm	120Å	10µm	C ₁₈	500g
30049023	300266	ZEOsphere 120 C ₁₈ / 10µm	120Å	10µm	C ₁₈	1kg
30049004	300292	ZEOsphere C ₈ 100Å / 10µm	100Å	10µm	C ₈	100g
30048994	302056	ZEOsphere C ₈ 100Å / 10µm	100Å	10µm	C ₈	500g
30049016	301153	ZEOsphere C ₈ 100Å / 10µm	100Å	10µm	C ₈	1kg
30049008	300293	ZEOsphere 120A C ₈ / 10µm	120Å	10µm	C ₈	100g
30048980	302059	ZEOsphere 120A C ₈ / 10µm	120Å	10µm	C ₈	500g
30048998	300277	ZEOsphere 120A C ₈ / 10µm	120Å	10µm	C ₈	1kg
30049028	302061	ZEOsphere 100 Diol / 10µm	100Å	10µm	Diol	100g
30048993	302062	ZEOsphere 100 Diol / 10µm	100Å	10µm	Diol	500g
30048979	302063	ZEOsphere 100 Diol / 10µm	100Å	10µm	Diol	1kg
30048977	302109	ZEOsphere 120 Diol / 10µm	120Å	10µm	Diol	100g
30048969	302110	ZEOsphere 120 Diol / 10µm	120Å	10µm	Diol	500g
30049018	302060	ZEOsphere 120 Diol / 10µm	120Å	10µm	Diol	1kg
30048981	302111	ZEOsphere 120 NH ₂ / 10µm	120Å	10µm	NH ₂	100g
30049032	302212	ZEOsphere 120 NH ₂ / 10µm	120Å	10µm	NH ₂	500g
30049011	302213	ZEOsphere 120 NH ₂ / 10µm	120Å	10µm	NH ₂	1kg
30048986	300866	ZEOsphere 300 C ₄ / 10µm	300Å	10µm	C ₄	100g
30048951	302068	ZEOsphere 300 C ₄ / 10µm	300Å	10µm	C ₄	500g
30049001	301284	ZEOsphere 300 C ₄ / 10µm	300Å	10µm	C ₄	1kg

Visit eu.fishersci.com for more information

Distributed by Fisher Scientific. Contact us today:

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk
 Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it
 Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl
 Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es
 Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

